



LONG-TERM CARE

# Total Hand Hygiene

CLEAN WITHOUT A DOUBT.

**ECOLAB**<sup>®</sup>

Everywhere It Matters.<sup>™</sup>

# Table of Contents

---

Nexa™ Concentrates

3

---

Nexa™ Dispensers

5

---

Nexa™ Hand Sanitizer Stations

8

---

Syncra™

9

---

Counter-Mount, Bulk & Tabletop


11

---

# Spread comfort, not infection.

In long-term care facilities, your priority is the health and safety of your residents. Implementing effective handwashing practices is the first line of defense against the spread of infection throughout your community. At Ecolab, we've developed a total hand hygiene system that combines effective chemistry with technology and training—making your hand cleaning habits efficient, safe† and simple.

## WHY FOCUS ON HANDWASHING?



### 1 TO 3 MILLION

serious infections occur every year in [long-term care] facilities.<sup>1</sup>

### 95% OF PEOPLE

don't wash their hands adequately.<sup>2</sup>

### THREE DAYS

is the amount of time some germs can live on moist surfaces.<sup>3</sup>

<sup>1</sup><http://www.cdc.gov/longtermcare/>

<sup>2</sup><http://www.flu.gov/planning-preparedness/hospitality/influenzaguidance.html>

<sup>3</sup><http://msutoday.msu.edu/news/2013/eww-only-5-percent-wash-hands-correctly/>

## AN EASIER WAY TO IMPROVE HAND HYGIENE.



Nexa™ Concentrates

Nexa™ Dispensers

Syncra™

Touch Point Service<sup>SM</sup>

**1 800 35 CLEAN // ECOLAB.COM**

To learn more about how Ecolab can help improve your compliance, safety and savings, contact your Ecolab Sales and Service Associate.

†Safe as long as products are used as directed.

# Nexa™ Concentrates


A sustainable solution  
that helps save money.

Ecolab's Nexa Concentrates provide sustainability and cost savings with uniquely designed refillable bottles for use with both Syncra™ and Nexa dispensers.



## Reduce waste.

Ecolab's sustainably designed system reduces packaging waste by up to 90 percent. The refillable design virtually eliminates product waste and makes refills a predictable part of your routine.<sup>1</sup>



## Save on shipping costs.

The concentrated formula of Nexa Concentrates mean you need to order less frequently, reducing shipping costs and saving on storage space. Fewer, smaller packages also helps reduce waste pickup.



## Wash with confidence.

Nexa Concentrates' refillable bottles feature a sealed system and in-chemistry preservatives to ensure safe<sup>2</sup> and effective application time-after-time.

# Nexa™ Concentrates



DISPENSER  
AND PRODUCTS



NEXA CONCENTRATES 1.3L SYSTEM	PART NUMBER	QUANTITY
1.3L Bottle Version for General Use and Sanitizer	9202-1202	1 Unit
1.3L Bottle Version for Antibacterial	9202-1204	1 Unit

REFILLABLE BOTTLES & EQUIPMENT	PART NUMBER	QUANTITY
1250 mL Refillable Bottle	9202-3330	4 Pack
	9202-3329	10 Pack
750 mL Refillable Bottle	9202-3328	4 Pack
	9202-3327	10 Pack
Nexa Bottle Caddy (Holds 6 Bottles)	9202-3356	1 Unit

LABELS FOR REFILLABLE BOTTLES	PART NUMBER	QUANTITY
Foam Hand Soap	9202-3317	10/Sheet
Antibacterial Foam Hand Soap	9202-3318	10/Sheet
Foam Hand Sanitizer	9202-3319	10/Sheet

PRODUCT	NAME	DESCRIPTION	PACK SIZE	PRODUCT CODE
General Use Hand Soaps	Concentrated Foam Hand Soap	Concentrated <b>Green Seal™</b> certified general use foam hand soap, that is ideal for public restrooms. It can be diluted up to 10:1 dilution rate. 	2-1.3 Liter	6100869
	Concentrated Mild Foam Hand Soap	Fragrance-free, concentrated <b>Green Seal™</b> certified general use foam hand soap, that is ideal for public restrooms. It can be diluted up to 10:1 dilution rate. 	2-1.3 Liter	6100875
Antibacterial Hand Soaps	Concentrated Antibacterial Foam Hand Soap	Concentrated FDA-regulated foaming antibacterial soap with a benzalkonium chloride active and a light citrus sent. It can be diluted up to a 5:1 dilution rate.	2-1.3 Liter	6100871
Hand Sanitizers	Concentrated Foam Hand Sanitizer	FDA food code-compliant and alcohol free concentrated foam hand sanitizer with a benzalkonium chloride active. It can be diluted up to a 10:1 dilution rate.	2-1.3 Liter	6100873



This product meets Green Seal™ Standard GS-41A based on effective performance and protective limits on VOCs and human & environmental toxicity. [GreenSeal.org](http://GreenSeal.org)



Concentrated Antibacterial Foam Hand Soap, Nexa Concentrates Dispenser

**1 800 35 CLEAN // ECOLAB.COM**

To learn more about how Nexa Concentrates can help you save money, contact your Ecolab Sales and Service Associate.

# Nexa™ Dispensers

Easy, efficient  
handwashing.

The Nexa Dispensing System provides efficiency, innovation and sustainability in one sleek, but adaptable, dispenser. Available in manual or touch-free.



## Use fewer batteries.

Energy saving touch-free dispenser requires fewer battery change-outs.<sup>1</sup>



## Inhibit bacterial growth.

Antimicrobial protection in the manual dispenser push bar continuously blocks the growth of bacteria, mold and mildew.<sup>2</sup>



## Optimize product use.

Improved HDPE<sup>3</sup> packaging is not only recyclable, it also helps to maximize product evacuation—dispensing up to 99 percent of product.



## Prevent costly waste.

Patented pump pulls excess product back in, reducing clogs and drips that result in costly product waste and preventing potential slips and falls.

<sup>1</sup>Over the course of 50 months at 50 activations per day, in comparative testing using Ecolab's Nexa Touch-Free dispenser and leading competitor's dispenser.

<sup>2</sup>Bacteria, mold and mildew that can cause degradation.

<sup>3</sup>HDPE: High-Density Polyethylene

DISPENSERS	BLACK	WHITE
1 Nexa Classic Touch-Free Dispenser (1250 mL)	9202-1192	9202-1193
2 Nexa Compact Touch-Free Dispenser (750 mL)	9202-1190	9202-1191
3 Nexa Classic Manual Dispenser (1250 mL)	9202-3091	9202-3093
4 Nexa Compact Manual Dispenser (750 mL)	9202-3087	9202-3088
Hand Sanitizer Station (Includes White Nexa Classic Touch-Free Dispenser) <b>SEE PAGE 8</b>	9202-6037	1 Unit

ACCESSORIES	PART NUMBER	QUANTITY
5 Nexa Security Enclosure (does not include Nexa dispenser)	9202-3448	1 Unit
Drip Tray (Black)	9202-2989	1 Unit
Drip Tray (White)	9202-2993	1 Unit
Lock and Key Kit	9202-3113	2 Keys/10 Locks
Individual Lock and Key	9202-3195	1 Key/1 Lock
6 Product Badges (Antimicrobial Hand Soap)	9202-3107	6 Badges
Product Badges (Hand Sanitizer)	9202-3106	6 Badges
Product Badges (Hair and Body Wash)	9202-3406	6 Badges
Product Badges (Hand Soap)	9202-3409	6 Badges
Product Badges (Moisturizer)	9202-3408	6 Badges

Continued on next page.



Touch-Free Dispensers

Manual Dispensers

Nexa Security Enclosure

Product Badges

\*Logo-only badge comes standard with each dispenser order.

**CUSTOMIZE YOUR OWN BADGE ONLINE:** [www.ecolab.com/handcaregraphicsbuilder](http://www.ecolab.com/handcaregraphicsbuilder)

**1 800 35 CLEAN // ECOLAB.COM**

To learn more about how Nexa dispensers can help increase efficiency and improve sustainability, contact your Ecolab Sales and Service Associate.



PRODUCT	NAME	DESCRIPTION	PACK SIZE	PRODUCT CODE	
General Use Hand Soaps	Foam Hand Soap	Green Seal™ certified gentle foaming hand soap that is ideal for public restrooms.		6-750 mL	6100708
				4-1250 mL	6100709
	Mild Foam Hand Soap	Green Seal™ certified dye-free, fragrance-free gentle foam hand soap ideal for public restrooms.		6-750 mL	6101194
				4-1250 mL	6101193
	Clean & Smooth	Green Seal™ certified lotion hand soap that provides rich lather to remove stubborn dirt while leaving hand fresh and soft.		6-750 mL	6100706
				4-1250 mL	6100707
Antibacterial Hand Soaps	Advanced Antibacterial Foam Hand Soap	Pre-lathered soap with a light floral scent and benzalkonium chloride active. Effective against most common bacteria that may cause illness in food service environments.		6-750 mL	6101090
				4-1250 mL	6101088
	Advanced Antibacterial Clean & Smooth	High-quality lotion soap that is gentle on hands with a benzalkonium chloride active for antimicrobial efficacy.		6-750 mL	6101421
				4-1250 mL	6101420
	Gentle Antibacterial Foam Hand Soap	Pre-lathered soap with a citrus scent and PCMX active. Effective against most common bacteria that may cause illness in food service environments.		6-750 mL	6100719
				4-1250 mL	6100720
	Medicated Foam Hand Soap	Antibacterial hand wash that meets FDA requirements for Healthcare Personnel Handwash. Contains aloe vera and emollients to help moisten and soothe the skin.		6-750 mL	6100717
				4-1250 mL	6100718
Hand Sanitizers	Foam Hand Sanitizer	Formulated, with 62% ethyl alcohol, to promote skin health and to reduce the spread of germs. 94% USDA certified biobased product.		6-750 mL	6100727
				4-1200 mL	6100728
	Foodservice Foam Hand Sanitizer	FDA food code-compliant (non-alcohol based). Contains skin conditioners and emollients to help moisturize and improve skin health. Ideal for use in the kitchen for foodservice settings.		6-750 mL	6100729
				4-1250 mL	6100730
	Gentle Hand Sanitizer Gel	Formulated with 70% w/w ethyl alcohol, killed 99.99% of pathogens within 15 seconds in in-vitro testing. Built with vitamin E, moisturizes and 5 emollients to help maintain skin moisture.		6-750 mL	6100725
				4-1200 mL	6100726
Instant Hand Sanitizer Gel	Food code compliant and more than 99.5% effective in in-vitro testing against a broad range of foodborne illness-causing bacteria. Contains a unique ingredient for superior after-feel.		6-750 mL	6100722	
			4-1200 mL	6100723	
Hand Lotions	Revitalizing Skin Lotion	Fragrance free, advanced formula designed to be non-greasy and quick absorbing providing great moisturization.		6-750 mL	6000062
Body Shampoo	Total Body Shampoo	A protein enriched, conditioning hair and body wash formulated to soothe hair and skin. Can be diluted for use as a no-rinse bed/bath or perineal cleanser.		6-750 mL	6000064
				4-1250 mL	6101363



This product meets Green Seal™ Standard GS-41A based on effective performance and protective limits on VOCs and human & environmental toxicity. [GreenSeal.org](http://GreenSeal.org)



These products are derived from renewable agricultural material. BioPreferred is a USDA program designed to increase the development and purchase of biobased products. For more information go to [www.biopreferred.gov](http://www.biopreferred.gov)

**1 800 35 CLEAN // ECOLAB.COM**

To learn more about how Nexa dispensers can help increase efficiency and improve sustainability, contact your Ecolab Sales and Service Associate.



# Nexa™ Hand Sanitizer Stations

Clean hands, wherever you need them.

Attractive and adaptable, Ecolab's Nexa hand sanitizer stations can be placed anywhere throughout your facility. Nexa waterless hand sanitizer delivers safe<sup>1</sup>, effective protection. The Nexa Touch-Free dispenser helps prevent the spread of germs. Uniquely designed product packaging helps maximize product evacuation—dispensing up to 99 percent of product.

## What you get when you order a Nexa hand sanitizer station:


<b>1</b>	<b>ONE NEXA CLASSIC TOUCH-FREE DISPENSER (IN WHITE)</b>
<b>2</b>	<b>ONE NEXA HAND SANITIZER STATION</b>
<b>3</b>	<b>TWO STANDARD GRAPHICS (OPTIONS SHOWN RIGHT)</b>



Hand Sanitizer Station (Includes White Nexa Classic Touch-Free Dispenser) 9202-6037

## DESIGN YOUR OWN HAND SANITIZER STATION GRAPHIC.

Make your hand sanitizer station stand out with your own logo and graphics! Visit the website below for design specifications and details about how to order.



CUSTOMIZE YOUR SANITIZER STATION AT [www.ecolab.com/handcaregraphicsbuilder](http://www.ecolab.com/handcaregraphicsbuilder)

**1 800 35 CLEAN // ECOLAB.COM**

To learn more, contact your Ecolab Sales and Service Associate.

# Syncra™

## Effective handwashing made simple.

A breakthrough innovation in hand hygiene, Syncra's automated touch-free water and soap delivery system provides more effective handwashing results. With a pre-programmed, three-step process for handwashing and rinsing, Syncra ensures consistent compliance with every wash.



### 3 EASY STEPS

- ✓ WATER & SOAP
- ✓ TIMED LATHER & SCRUB
- ✓ RINSE



### Improve handwashing compliance.

Syncra combines easy-to-use technology with safe<sup>1</sup>, effective germ-killing power to help increase handwashing compliance.



### Eliminate the "ick factor."

With touch-free activation, Syncra bypasses germs that can grow on traditional faucet handles and soap dispensers.




### Reduce water use.

Syncra's timed delivery saves on excessive water use—dispensing just the right amount of water for a thorough clean. This can reduce your water usage by up to 50 percent<sup>2</sup>.

<sup>1</sup>Safe as long as products are used as directed.

<sup>2</sup>Water reduction calculated vs. a traditional faucet which runs based on FDA Food Code recommended wash time.

EQUIPMENT	PART NUMBER	QUANTITY
4" Wall Mount Faucet with Dispenser Enclosure	9202-1225	1 Unit
Center Deck Mount Faucet with Dispenser Enclosure	9202-1226	1 Unit
4" Deck Mount Faucet with Dispenser Enclosure	9202-1227	1 Unit
8" Deck Mount Faucet with Dispenser Enclosure	9202-1228	1 Unit

PRODUCT	NAME	DESCRIPTION	PACK SIZE	PRODUCT CODE
General Use Hand Soaps	Foam Hand Soap	Green Seal™ certified gentle foaming hand soap that is ideal for public restrooms. 	4-1250 mL	6100865
Antibacterial Hand Soaps	Advanced Antibacterial Foam Hand Soap	Antibacterial foam soap with a benzalkonium chloride active and a light floral scent.	4-1250 mL	6101087



This product meets Green Seal™ Standard GS-41A based on effective performance and protective limits on VOCs and human & environmental toxicity. [GreenSeal.org](http://GreenSeal.org)



\*Syncra Faucet

**1 800 35 CLEAN // ECOLAB.COM**

To learn more about how Syncra can help improve your compliance, safety and savings, contact your Ecolab Sales and Service Associate.

\*Sink not included with Syncra Faucet unit

# Counter-Mount, Bulk & Tabletop

An array of versatile hand hygiene options.

- 1 COUNTER-MOUNT
- 2 TABLETOP
- 3 BULK



Bulk Dispenser, EQUI-MILD™ Foam, Counter-Mount Dispenser




Designed for better handwashing.

Help prevent the spread of germs with Ecolab's counter-mount, touch-free dispenser. This easy-to-use system provides reliable, low-maintenance dispensing of Ecolab foam hand soap products. Designed to meet your aesthetic needs.



Clean hands anywhere, anytime.

For flexible hand hygiene where and when it's needed, choose Ecolab tabletop solutions. These conveniently sized products are packaged for use without dispensers so they can be placed wherever necessary. For reliable sanitizing on the go, Ecolab offers a variety of hand hygiene products packaged for personal use.



Flexibility with bulk solutions.

Ecolab offers a variety of hand hygiene products in bulk quantities for use in refillable dispensers.

# Counter-Mount

DISPENSER  
AND PRODUCTS



DISPENSERS	PART NUMBER	QUANTITY
Counter Mount Dispenser	9263-1012	1 Unit

PRODUCT	NAME	DESCRIPTION	PACK SIZE	PRODUCT CODE
General Use Hand Soaps	Foam Hand Soap	A foaming hand soap that is ideal for public restrooms. It is enriched with emollient and moisturizers to keep hands smooth and is gentle enough for repeat use.	2-1600 mL	6100122
Antibacterial Hand Soaps	Advanced Antibacterial Foam Hand Soap	Antibacterial foam soap with a benzalkonium chloride active and a light floral scent.	2-1600 mL	6101086

Continued on next page.



Counter Mount Dispenser, Digiclean Mild Foam Hand Soap

**1 800 35 CLEAN // ECOLAB.COM**

To learn more, contact your Ecolab Sales and Service Associate.


PRODUCT	NAME	DESCRIPTION	PACK SIZE	PRODUCT CODE
General Use Hand Soaps	Endure Sensitive Skin and Hand Soap	Dye-free, fragrance free formula is enriched with skin conditioners and emollients to help protect against irritation and keep skin soft and smooth.	12-500 mL	6000116
	Antibacterial Hand Soaps	EQUI-MILD™ Foam Antibacterial Hand Soap	A gentle foam soap with a light floral scent. The Triclosan-free formula uses benzalkonium chloride (quat) active to be effective against most common bacteria that may cause illness in food service environments.	12-500 mL
		Medi-Stat	A mild Healthcare Personnel Hand Wash containing 0.5% Chloroxylenol (PCMX) as an antimicrobial active. Contains aloe vera to help moisten and soothe skin during frequent use.	12-540 mL
Hand Sanitizers	Quik-Care Foam Hand Sanitizer	Dye-free, fragrance-free formula provides fast antimicrobial efficacy against a wide range of microorganisms with reduced risk of irritation and allergic reactions.	12-500 mL	6032105
	Advanced Gel Hand Sanitizer	A waterless antimicrobial hand rinse with an emollient-rich lotion formula to help moisturize and repair skin. Gives the added protection by helping control the spread of nosocomial infections when conditions prevent routine hand washing. Advanced Gel Hand Sanitizer is compatible with CHG, latex gloves and lotions.	12-540 mL	6030370
			72-4 oz	6030360
			144-1.25 oz	6030392
	Gel Hand Sanitizer	Food-code compliant hand sanitizer gel that reduces the level of bacteria on hands. For use between soap and water handwashing. This formula is more than 99.5% effective in in-vitro testing against a broad range of foodborne illness-causing bacteria.	12-540 mL	6000004
			72-4 oz	6000003
144-1.25 oz			6030352	
Hand Lotion	Revitalizing Skin Lotion	Fragrance free, advanced formula designed to be non-greasy and quick absorbing providing great moisturization.	72-4 oz	6059384
			12-540 mL	6059323
Body Shampoo	Total Body Shampoo	A protein enriched, conditioning hair and body wash formulated to soothe hair and skin. Can be diluted for use as a no-rinse bed/bath cleanser.	72-4 oz	6021162
			12-1000 mL	6067132



Gel Hand Sanitizer, EQUI-MILD™ Foam Antibacterial Hand Soap

DISPENSERS	PART NUMBER	QUANTITY
Manual Bulk Refillable Dispenser	9263-1970	1 Unit

ACCESSORIES	PART NUMBER	QUANTITY
Manual Bulk Dispenser Drip Tray	9202-1086	12 Pack

PRODUCT	NAME	DESCRIPTION	PACK SIZE	PRODUCT CODE
General Use Hand Soaps	Clean & Smooth	Green Seal™ certified lotion soap that provides rich, moisturizing lather to remove stubborn dirt while leaving hands fresh and soft. 	4-1 gal	6115875
	Liquid Hand Soap	Lotion soap with light, fresh fragrance. Extra mild formula is ideal for frequent handwashing.	4-1 gal	6100313
Antibacterial Hand Soaps	Advanced Antibacterial Clean & Smooth	Is an economical, high-quality antibacterial lotion soap that is gentle on hands and effective against germs. It is designed to keep employees hands clean and sanitary and to reduce the threat of foodborne illness.	4-1 gal	6101419
Body Shampoo	Total Body Shampoo	A protein enriched, conditioning hair and body wash formulated to soothe hair and skin. Can be diluted for use as a no-rinse bed/bath or perineal cleanser.	4-1 gal	6100628



This product meets Green Seal™ Standard GS-41A based on effective performance and protective limits on VOCs and human & environmental toxicity. [GreenSeal.org](http://GreenSeal.org)



# DigiGlo

## HAND WASH TRAINING

Manual Bulk Refillable Dispenser, Advanced Antibacterial Clean & Smooth

PRODUCT	NAME	DESCRIPTION	PACK SIZE	PRODUCT CODE
Training Kit	DigiGlo Fluorescent Indicator Demo Kit	Helps teach proper handwashing techniques. Each kit contains two 4oz bottles of the fluorescent indicator lotion (~100 uses), a black light flashlight with batteries, and an instruction sheet, all in a plastic case.	1 Unit	6100903



DigiGlo Fluorescent Indicator Demo Kit

# Touch Point Service<sup>SM</sup>

Expert guidance with a personal touch.

When it comes to training and service, we've got you covered with Ecolab's Total Hand Hygiene System and Touch Point Service<sup>SM</sup>.



## WE'LL INSTALL IT FOR YOU... FOR FREE.

Complete installation service managed by Ecolab experts

Assist locations in following fire code specifications and ADA guidelines



## EXPERT TRAINING FOR EFFECTIVE USE.

Training on proper product and equipment usage and maintenance

Visual aids provide ongoing education and reminders of proper procedures



## SOLUTIONS DESIGNED FOR YOU.

Customized product selection, tailored to your needs

Consistent implementation across locations

Advice on product placement in each location

Knowledgeable suggestions on product level management

Ongoing inventory monitoring



## ENGAGE YOUR STAFF WITH CONTINUING EDUCATION.

Educational materials that help inform and encourage staff to keep hand hygiene top of mind

Posters and point-of-use materials that engage employees and guests



## WE'LL BE THERE TO HELP YOU KEEP IMPROVING.

Regular service visits to evaluate and improve each location's hand hygiene program

Ongoing consultations to boost results and efficiency



To order,  
contact your  
Ecolab Sales and  
Service Associate.